

View Online at <https://aerobasegroup.com/nsn/6630-01-469-3390>

**Special Features:**

Quality control solution can be introduced into analyzer by syringe, capillary, or directly from vial; solution tests analyzer for low limits of ph (hydrogen ion concentration), po<sub>2</sub> (partial pressure oxygen), and pco<sub>2</sub> (partial pressure carbon dioxide); if analyzer fails test, patient blood sample will be refused; packaging, 30 - 1.7 milliliter ampoules

**Functional Description:**

Used to calibrate opti-1 portable blood gas analyzer

**Type/model Designation:**

Opti-trol level 1

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No