

View Online at <https://aerobasegroup.com/nsn/6630-01-469-3390>

Special Features:

Quality control solution can be introduced into analyzer by syringe, capillary, or directly from vial; solution tests analyzer for low limits of ph (hydrogen ion concentration), po₂ (partial pressure oxygen), and pco₂ (partial pressure carbon dioxide); if analyzer fails test, patient blood sample will be refused; packaging, 30 - 1.7 milliliter ampoules

Functional Description:

Used to calibrate opti-1 portable blood gas analyzer

Type/model Designation:

Opti-trol level 1

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No