

View Online at <https://aerobasegroup.com/nsn/6640-01-347-1092>

**Product Name:**

Priming solution automated analysis 16 ml bottle 5's

**General Description:**

116 ml bottle; prepares the sample line before aspirating coulter's latron control in order to ensure a background free, electrostatic free pathway for the control; used with coulter's stks, maxm and vcs clinical hematology instruments; unit qty: 5 per pkg; this item is usually shipped with 1 year shelf life, expires 30 days after opening; latron primer prepares the tubing and instrument components for the latron control for coulter hmx hematology analyzer; refrigeration is required at 2-8 deg.C

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No