

View Online at <https://aerobasegroup.com/nsn/6640-01-468-3542>

**End Application:**

Used to make internal structures visible under magnification by clearing out soft tissue

**Product Name:**

Clearing fluid specimen lacto-phenol sol w/glacial acetic acid

**Special Features:**

Clearing fluid, specimen lacto-phenol solution w/glacial acetic acid. Specimens may remain in fluid for prolonged periods of time or may be heated to speed the clearing process

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No