

View Online at <https://aerobasegroup.com/nsn/6640-01-509-0794>

**Product Name:**

Decontaminant surface dna-away 250ml bt

**Special Features:**

Specially formulated surface decontaminant eliminates unwanted dna and dnase from glassware and plasticware without affecting subsequent cna samples; use before performing pcr or working with dna for a dna-free work area; degrades dna more effectively than autoclaving; plus it's quicker; simply apply the ready-to use formula to the surface or labware to be decontaminated; wipe dry or rinse clean; suitable for gel boxes, pipettors, benchtops, thermal cyclers or any other apparatus used when working with dna; 8.5oz. (250ml) bottle

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No