

View Online at <https://aerobasegroup.com/nsn/6640-01-564-0276>

**Furnished Items:**

Unisex collection swabs; transport media

**End Application:**

6550-01-522-4071, test kit, chlamydia and neisseria detection

**Product Name:**

Specimen collection kt chlamydia gc endocervial male urethra 50s

**Special Features:**

Kit consists of collection swabs and transport media for collection of endocervical and male urethra specimens for the detection of chlamydia trachomatis and neisseria gonorrhoea via pcr testing

**Supplementary Features:**

After collection, specimen must be transported to the laboratory at 2 deg. To 30 deg.C if they will be tested within 60 days of collection; specimen can be frozen at neg. 20 degrees to neg. 70 degrees c to allow testing up to 90 days

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A23800