

View Online at <https://aerobasegroup.com/nsn/6840-01-485-1910>

Special Features:

Cidex opa solution test strips are semi-quantitative chemical indicators for use in determining the concentration of ortho-phthalaldehyde, the active ingredient in cidex opa; the test strip is composed of paper impregnated with sodium sulfite and ph-sensitive dye and buffers; 2 year shelf life; 4 packs of 250, 1000 strips per case; do not freeze

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No