

View Online at <https://aerobasegroup.com/nsn/7030-01-600-9011>

**End Application:**

Bd cytometric bead array (cba) analysis

**Product Name:**

Fcap array v3.0 software

**Special Features:**

Intuitive layout of samples and assignment to data files. Standard dilution calculator for fast dataentry. Debris filtering and manual clustering options. Results viewing in tabular or graphical format by sample or by analyte. Campatiblility with fcs 2.0 and fcs 3.0 data files. Automated workflow for the bd facsverse flow cytometer. Contains installation dvd, hardware key, user guide.

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A23800