## NSN 4110-01-157-5898

Refrigerated Bath Circulators - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/4110-01-157-5898

## **General Description:**

1bath uses a fully hermetic compressor which is sealed permanetly, operating at temp of +40 degrees c; a control panel switch allows the
compressor to be shut off on all baths for fast heat-up; z pump can be circulated through two closed loop system

1bath uses a fully hermetic compressor which is sealed permanetly, operating at temp of +40 degrees c; a control panel switch allows the
compressor to be shut off on all baths for fast heat-up; z pump can be circulated through two closed loop system
Shelf Life:
N/a
Unit Of Measure:
Demilitarization:
No