Coolanol-slip Ring Seal Bladder - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/5342-01-317-6862

End Application:

Lantirn targeting pod

General Description:

Material-bag and inlet tube flourosilicone per mil-r-25988, type 2, cl 1 gr 60 or ams 33260, swelling in coolanol 25r shall not exceed 8% when tested per astm d471, volume resisitivity per astm d991 shall be less than 5000 ohm cm; tube hole od 1.250 in., id 0.870 to 0.900 in. W/0.050 in. Porm 0.005in. Wall thickness within 1.250 in. Diameter & 0.050 in. All other places; reservoir to hold a of 200 cc of coolanol when constrained

Shelf Life:

N/a

Unit Of Measure:

Demilitarization:

Yes - demil/mli