Impregnated Cast Bandage - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6510-01-397-6915

## Setting Time:

Between 3 minutes and 4 minutes

## **Features Provided:**

Perforated

## **Product Name:**

Orthoplast ii bandage cast impregnated molds and contours 2 sheets

## **Special Features:**

Transpolyisoprene thermoplastic splinting material; use warm water or dry heat to form splint; molds and contours w/o wrinkling; material will bond to itself w/additional heat and it's edges are self-finishing; smooth; 24 in. X 36 in.; primarily designed for upper extremities; wax coated surface prevents accidental self-adhesion; unit qty is 2 sheets per package

Shelf Life:

N/a

---

Unit Of Measure:

Demilitarization:

No

Fiig:

т00600