

View Online at <https://aerobasegroup.com/nsn/6515-00-656-1021>

Current Type:

Ac

Case Material:

Plastic

Accessory Components And Quantity:

Adapter, intradermal nozzle cleaning, 1; adapter, needle, 3; adapter, subcutaneous nozzle cleaning, 1; ball, nylon, 3/32 in. Dia, 8; ball, nylon, 3/16 in. Dia, 3; ball, cres, 3/16 in. Dia, 1; cap, inj. Nozzle, 1; cloth, polishing, 1; nozzle, injection, subcutaneous, 1; o-rings, assorted sizes, 6; o-ring, polytetrafluoroethylene, size 5, 3; quad seals, assorted sizes, 8; ring, packaging, back-up, 2; spring, firing, cocking and check valve, 3; spring, main, 1; spring, outlet valve, 3; tool, maintenance, 1; tube, vaccine feed and vent, 3

Case:

Provided

Control Type:

Trigger and trigger

Electric Motor Horsepower Rating:

0.083

Dosage Range:

+0.1 to 1.0 cubic centimeters

Nozzle Tip Material:

Steel, corrosion resisting

Tip Orifice Insert Material:

Sapphire

Control Function:

Loading and injection

Pump Operating Pressure:

1200.0 pounds per square inch

Quick Disconnect Hydraulic Hose Coupling Connection:

Included

Bottle Holder:

Included

Product Name:

Hypodermic injection apparatus jet automatic 115v 50/60 hz ac

Exterior Container Weight:

35.00 pounds

Voltage In Volts:

115.0

Exterior Container Cubic Measure:

1.00 cubic feet

Exterior Container Quantity:

1

Special Features:

FORR-AP-IDL-YINJ

Frequency In Hertz:

50.0 or 60.0

Unit Of Issue Cube:

1.200 cubic feet

Unit Of Issue Weight:

27.00 pounds

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T107-c