

View Online at <https://aerobasegroup.com/nsn/6515-00-910-0097>

**Case Material:**

Plastic

**Accessory Components And Quantity:**

Adapter, intradermal nozzle cleaning, 1; adapter, needle, nylon, 3; adapter, subcutaneous nozzle cleaning, 1; ball, nylon, 3/32 in. Dia, 8; ball, nylon, 3/16 in. Dia, 3; cap, injection nozzle, 1; cloth, polishing, 1; needle, vaccine feed, 3; o-ring, sizes 1, 5, 8, 2 ea; o-ring, polytetrafluoroethylene, size 5, 3; ring, back-up, polytetrafluoroethylene, 2; seal, quad, sizes 8, 15, 2 ea; spring, firing and check valve, 3; spring, outlet valve, 3; tube, feed and vent, 3

**Case:**

Provided

**Control Type:**

Trigger and trigger

**Err Default:**

Err default

**Err Default:**

Err default

**Err Default:**

Err default

**Err Default:**

Err default and err default

**Err Default:**

Err default

**Err Default:**

Err default

**Err Default:**

Err default

**Product Name:**

Hypodermic injection apparatus jet automatic: foot operated

**Special Features:**

For rapidly injecting subcutaneously and/or intramuscularly metered doses; medicinal head shall be sterilizable and have a non-recessed tip

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T107-c