

View Online at <https://aerobasegroup.com/nsn/6515-01-332-8950>

Nonsupply Items And Quantities:

Consists of arterial line 148 inches, and venous line 109 inches; 8mm i.D. Pump segment, luer lock connectorsthroughout, clamps, negative pressure pillow, saline administration port, priming ports with caps, heparin infusion line, post-pump arterial-pressure monitoring chamber, venous bubble trap, and injection sites. This set is specially designed for high-efficiency hemodialysis where blood flow rates and pressures may be higher

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A23800