

View Online at <https://aerobasegroup.com/nsn/6515-01-332-8950>

**Nonsupply Items And Quantities:**

Consists of arterial line 148 inches, and venous line 109 inches; 8mm i.D. Pump segment, luer lock connectorsthroughout, clamps, negative pressure pillow, saline administration port, priming ports with caps, heparin infusion line, post-pump arterial-pressure monitoring chamber, venous bubble trap, and injection sites. This set is specially designed for high-efficiency hemodialysis where blood flow rates and pressures may be higher

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

A23800