

View Online at <https://aerobasegroup.com/nsn/6515-01-397-4056>

End Application:

COMP-AT-IBL-EWIT

Unit Package Quantity:

1

Unit Package Type:

Any acceptable

Component And Quantity:

C/o trocar, 1; zero adjustment tool, 1; catheter with pressure transducer and connector, 1

Features Provided:

Disposable and sterile

Product Name:

Monitoring kit post-craniotomy subdural pressure disp sterile10s

Special Features:

For tunneling the catheter beneath the scalp; kit used to directly monitor post-operative subdural pressure during emergency or elective surgery

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T107-a