

View Online at <https://aerobasegroup.com/nsn/6515-01-407-6583>

**Unit Package Quantity:**

1

**Sterility:**

Nonsterile

**Disposition After Initial Use:**

Reusable

**Special Features:**

Shoulder type, compression; consists of 3 parts: cuff, collar and tube; cuff will cover the shoulder with pressurized ice water, cooler will hold enough water and ice for 6 to 8 hours, tube will exchange the water between cooler and cuff; to minimize hemarthrosis and swelling, and cold to minimize pain

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

A23800