

View Online at <https://aerobasegroup.com/nsn/6515-01-407-6583>

Unit Package Quantity:

1

Sterility:

Nonsterile

Disposition After Initial Use:

Reusable

Special Features:

Shoulder type, compression; consists of 3 parts: cuff, collar and tube; cuff will cover the shoulder with pressurized ice water, cooler will hold enough water and ice for 6 to 8 hours, tube will exchange the water between cooler and cuff; to minimize hemarthrosis and swelling, and cold to minimize pain

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A23800