## NSN 6515-01-442-7810

T133-f

Patient Vital Signs Monitor - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6515-01-442-7810

| View Online at https://aerobasegroup.com/nsn/6515-01-442-7810   |
|---|
| Current Type:   |
| Ac/dc   |
| Product Name:   |
| Monitor patient vital signs 110/220v 50/60hz ac/dc  |
| Voltage In Volts:   |
| 110.0 and 220.0 and 12.0  |
| Special Features:   |
| C/o: 90309 pc scout tranportable/bedside monitor; 3 trace el display; dual channel internal recorder; 90351-1d expansion housing for the      |
| 90309 scout to provide spo2 and adult/neonatal nibp; 90470-07-21 multi-lead adult/neonatal ecg module with esis; two invasive blood           |
| pressures; two temperatures; multiview ii arrythmia review/edit and st segment analysis; 012-0291-01 pcms accessory kit; adult and            |
| neonatal nibp hose; one box of neonatal tru-cuff disposable nibp cuffs; one box of adult tru-cuff disposable nibp cuffs (please specify size) |
| one battery p/n 146-0018-00; one ac converter (must specify 110 volt p/n 650-0379-00 or 220 volt p/n 650-0391-00).                            |
| Frequency In Hertz:   |
| 50.0 and 60.0   |
| Supplementary Features:   |
| Note to requisitioner: the p/n 999-nsn2-00 was assigned by spacelabs to represent a configured system requested by the military services      |
| The component listing for the system can be found in "features".  |
| Shelf Life:   |
| N/a   |
| Unit Of Measure:  |
| <del></del>   |
| Demilitarization:   |
| No  |
| Fiig:   |