

View Online at <https://aerobasegroup.com/nsn/6515-01-477-4191>

**General Description:**

The ox174 dna/hae iii markers are made by digesting doubled-stranded 0x174 dna to completion with hae iii, the 11 dna fragments range from 72 and 1, 353 base pair pairs, blue/orange loading dye (6x) is provided, size 50ug, concentration 1ug/ul, recommended loading 1ul, typical number of lanes 50, range (bp) 72-1, 353, number of bands 11, storage condition minus 20 degrees c, storage buffer 10mm tris-hcl, ph 7.5, 1 mm edta

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No