

View Online at <https://aerobasegroup.com/nsn/6515-01-486-0188>

**Current Type:**

Ac

**Frequency Rating:**

Between 50.0 hertz and 60.0 hertz

**End Application:**

Used in assemblages where electrosurgical support is required

**Installation Design:**

Mobile

**Voltage Rating:**

Between 115.0 volts and 230.0 volts

**Unit Type:**

Solid state

**Mounting Type For Which Designed:**

Table

**Blended Output Current Feature:**

Included

**Product Name:**

Electrosurgical apparatus mobile table mnt

**Phase:**

Single

**Special Features:**

Weights less than 15 lbs; includes a return electrode quality monitoring system (recqms) that continually monitors impedance levels and deactivates the generator if a fault in the patient/return electrode contact is detected

**Supplementary Features:**

Surgistat ii electrosurgical generator or electrosurgical apparatus, compact

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T107-e