NSN 6515-01-566-4033

Hand Operated Resuscitator - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6515-01-566-4033

Adult resuscitation Accessory Components And Quantity: Silicone bag 1600ml 1; pt valve w/o relief valve 1; mask 3/4 w/mask cover 1; intake/reservoir valve 1; reservoir bag 2600 1 Disposition After Initial Use: Reusable Material And Location: Silicone bag Operation Method: Hand pump Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No Fig:	End Application:
Silicone bag 1600ml 1; pt valve w/o relief valve 1; mask 3/4 w/mask cover 1; intake/reservoir valve 1; reservoir bag 2600 1 Disposition After Initial Use: Reusable Material And Location: Silicone bag Operation Method: Hand pump Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Adult resuscitation
Disposition After Initial Use: Reusable Material And Location: Silicone bag Operation Method: Hand pump Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Accessory Components And Quantity:
Reusable Material And Location: Silicone bag Operation Method: Hand pump Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Silicone bag 1600ml 1; pt valve w/o relief valve 1; mask 3/4 w/mask cover 1; intake/reservoir valve 1; reservoir bag 2600 1
Material And Location: Silicone bag Operation Method: Hand pump Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Disposition After Initial Use:
Silicone bag Operation Method: Hand pump Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Reusable
Operation Method: Hand pump Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Material And Location:
Hand pump Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Silicone bag
Product Name: Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Operation Method:
Reusable resuscitator Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Hand pump
Special Features: Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Product Name:
Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive Shelf Life: N/a Unit Of Measure: Demilitarization: No	Reusable resuscitator
Shelf Life: N/a Unit Of Measure: Demilitarization: No	Special Features:
N/a Unit Of Measure: Demilitarization: No	Standard adult over 55 pounds; compliance with essential requirements of council directive 93/42 eec; medical device directive
Unit Of Measure: Demilitarization: No	Shelf Life:
Demilitarization: No	N/a
No	Unit Of Measure:
No	
	Demilitarization:
Fiig:	No
	Fiig:
T133-k	T133-k