

View Online at <https://aerobasegroup.com/nsn/6515-01-567-5587>

General Description:

Ready to use density gradient medium, 2 concentrations; 50 pct upper layer/90 pct lower layer; 12 by 6ml vial kit, 18 determinations; 2 by 100ml vial kit, 50 determinations; 2 yr shelf life; store at 2 to 8 degrees c; separates motile sperm from seminal fluid by density centrifugation; removes contaminants; sterile; herpes buffered; not for co2 incubator unless tightly capped; no dilution required, ready to use

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No