

View Online at <https://aerobasegroup.com/nsn/6515-01-576-6878>

Special Features:

Bead-based immunoassay used to qualitatively detect igg antibodies to hla class i and ii antigens; consists of 120ul hli i beads, 170ul conjugate concentrate, 30ml wash buffer, 80ul positive and negative control serum all har mful if swallowed; all components are liquid; hla beads in a pbs-based buffer are light blue, conjugate concetrate is pink, wash buffer is clear and slightly yellow, positive and negative control serum is clear and slightly yellow to amber

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No