

View Online at <https://aerobasegroup.com/nsn/6515-01-576-6884>

**Special Features:**

Bead-based immunoassay used to qualitatively detect igg antibodies to hla class ii antigens; consists of 120ul hli i beads, 170ul conjugate concentrate, 30ml wash buffer, 80ul positive and negative control serums all harmful if swallowed; all components are liquid; hla beads in a pbs-based buffer are light blue; conjugate concentrate is pink, wash buffer is clear to slightly yellow and positive/negative serum is clear and slightly yellow to amber

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No