## NSN 6515-01-580-6551

Surgical Bone Cement - Page 1 of 1

No Fiig: A535p0



View Online at https://aerobasegroup.com/nsn/6515-01-580-6551 **End Application:** Orthobiologies and surgery **Physical Form:** Paste **Quantity Within Each Unit Package:** 1.0 cubic centimeters Specific Usage Design: Bone graft substitute **Features Provided:** Separate catalyst **Product Name:** Osteoinductive bone graft substitute; dynagraft ii **Special Features:** Ddemineralized bone matric in reverse phase medium-donor bone marrow aspirate; type i college (bovine achilles tendon); beta-tcp Material: Cement **Special Test Features:** Federal (usa) law restricts this device to sale by or on the order of a physician **Supplementary Features:** Do not use after opening; package sterilized with radiation; do not refrigerate; store at 10 to 30 degrees cel; comes as putty-thickens at body temperature; mixed with donor tissue at time of grafting Shelf Life: N/a **Unit Of Measure:** Demilitarization: