NSN 6515-01-598-6630

Orthopedic External Fixation Set - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6515-01-598-6630

Component And Quantity:

Delta rod to rod coupling dia5/8/11 rods and 5mm diameter pins, 16; multiplanar rod to rod coupling, 6; rod coupler 30 degrees angled, 4; apex half pins 5mmx180mm 50mm 5mm thread, 26; apex half pins 5mmx150mm 40mm 5mm thread, 12; apex half pins 5mmx150mm 40mm 4mm thread, 12; apex half pins 5mmx120mm 27mm 3mm thread, 25; apex transfixing pins 5mmx300mm 40mm, 2; removeable thumbwheel, 2; 5 hole multi-pin clamp, 16; 11mm diameter posts 30 degrees angled, 16; soft tissue protector, 2; trocar, 2; 7mm spanner wrench, 2; 7mm t-wrench dia4/5/6mm pin driver, 1; stabilization/reduction wrench, 2; ao/jacobs quick release apex 4/5/6mm chuck, 2; apex drill brace 3/4/5/6mm, 1; carbon connecting rod 8mm diax250mm, 2; carbon connecting rod 8mm diax400mm, 4; carbon connecting rod 11mm diameter x300mm, 2; carbon connecting rod 11mm diameter x400mm, 4; semi-circular rod 11mm diax400mm, 2

Features Provided:

Sterile

Product Name:

External fixation tub system level 3/4/5 hoffman iii

Special Features:

External fixation tub system, level 3/4/5; sterile; consists of interchangeable components and instruments; snap-fit mechanism facilitates assembly and temporary fixation during reduction maneuvers; clamps, posts, and couplings are made from a combination of aluminum and austentic steel; carbon connecting rods are coated with electrically insullated vectran; 21.6 in l x 10.2 in w x 4.3 in d; approx. Weight: 20.1 lb; consists of three levels: base, lower insert, and upper insert

_					
S	hΔ	11		ife	
J			_	116	

N/a

Unit Of Measure:

Demilitarization:

No

Fiig:

T107-a