NSN 6515-01-619-7327

Catheter Connector - Page 1 of 1

No



View Online at https://aerobasegroup.com/nsn/6515-01-619-7327

:								
nector 100s	•							
es:								
heter, neut	ral displacement;	sterile;	flow rate:	165 ml/min;	75 psig back pressure	; priming vol:	0.05 ml;	box of 100
re:								
n:								
r	nector 100s es: neter, neut	nector 100s es: neter, neutral displacement; re:	nector 100s es: neter, neutral displacement; sterile; re:	nector 100s es: neter, neutral displacement; sterile; flow rate: re:	nector 100s es: neter, neutral displacement; sterile; flow rate: 165 ml/min; re:	nector 100s es: neter, neutral displacement; sterile; flow rate: 165 ml/min; 75 psig back pressure re:	nector 100s es: neter, neutral displacement; sterile; flow rate: 165 ml/min; 75 psig back pressure; priming vol: ee:	nector 100s es: neter, neutral displacement; sterile; flow rate: 165 ml/min; 75 psig back pressure; priming vol: 0.05 ml; re: