

View Online at <https://aerobasegroup.com/nsn/6520-00-044-1936>

**Unit Package Quantity:**

96

**Physical Form:**

Gel

**Unit Package Type:**

Package

**Setting Time:**

Between 8 minutes and 10 minutes

**Workable Consistency Temperature Rating:**

Between 109.0 degrees fahrenheit and 122.0 degrees fahrenheit

**Hydrocolloid Base Type:**

Agar

**Product Name:**

Impression material dental cartridge form hydrocolloid agar 48s

**Exterior Container Weight:**

24.00 pounds

**Exterior Container Cubic Measure:**

0.67 cubic feet

**Exterior Container Quantity:**

24

**Special Features:**

C/o 48 glass or plastic cartridges filled w/hydrocolloid gel and forty eight 19 or 21 gage disposable dispensing needles

**Intermediate Package Quantity:**

12

**Unit Of Issue Cube:**

0.028 cubic feet

**Supplementary Features:**

FORU-SE-W/C-ARTR

**Unit Of Issue Weight:**

1.00 pounds

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T243-d