NSN 6520-00-044-1936

Dental Impression Material - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6520-00-044-1936 **Unit Package Quantity: Physical Form:** Gel **Unit Package Type:** Package **Setting Time:** Between 8 minutes and 10 minutes **Workable Consistency Tempurature Rating:** Between 109.0 degrees fahrenheit and 122.0 degrees fahrenheit **Hydrocolloid Base Type:** Agar **Product Name:** Impression material dental cartridge form hydrocolloid agar 48s **Exterior Container Weight:** 24.00 pounds **Exterior Container Cubic Measure:** 0.67 cubic feet **Exterior Container Quantity:** 24 **Special Features:** C/o 48 glass or plastic cartridges filled w/hydrocolloid gel and forty eight 19 or 21 gage disposable dispensing needles **Intermediate Package Quantity:** 12 **Unit Of Issue Cube:** 0.028 cubic feet **Supplementary Features:** FORU-SE-W/C-ARTR **Unit Of Issue Weight:** 1.00 pounds Shelf Life: N/a **Unit Of Measure: Demilitarization:**

No

Fiig:

T243-d