

View Online at <https://aerobasegroup.com/nsn/6520-00-702-5557>

Physical Form:

Powder and liquid

Shade Source:

Ld caulk co

Shade Identification:

Light

Furnished Items And Quantity:

21g powders, 50; 430ml liquid, 1; 120ml liquid, 1; 10g gingival tone powders, 2; accessories

Usage Design:

Denture impression tray fabrication

Weight:

21.0 grams

Shade Name:

Characterized lucitone

Curing Process Type:

Heat

Special Features:

Lucitone characterized; denture base resin (c-25)

Color:

Pink, light

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d