Dental Impression Material - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6520-01-020-1428

Physical Form:
Stick
Setting Time:
Between 50 seconds and 80 seconds
Weight:
0.250 pounds
Physical Form Color:
Red
Workable Consistency Tempurature Rating:
Between 53.3 degrees celsius and 56.7 degrees celsius
Product Name:
Impression material dental sticks red 50/80 seconds set 0.25 lb
Specification Or Standard:
I type and 2 form
Departure From Cited Document:
As modified by aguc
Shelf Life:
N/a
Unit Of Measure:
Demilitarization:
No
Fiig:
T243-d