

View Online at <https://aerobasegroup.com/nsn/6520-01-272-6067>

**Length:**

10.000 feet

**Physical Form:**

Cord

**Disposition After Initial Use:**

Disposable

**Specific Use:**

Knitted cord used to retract the gingival area before placing impression material

**Product Name:**

Gingival retraction cord non-impregnated dental sz0 10' cotton

**Special Features:**

Size 0; used w/astringent liquid, pn 0111; recommended for anterior teeth; used as bottom cord where two cord technique is desired

**Material:**

Cotton

**Material:**

Cotton

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T243-e