

View Online at <https://aerobasegroup.com/nsn/6520-01-272-6067>

Length:

10.000 feet

Physical Form:

Cord

Disposition After Initial Use:

Disposable

Specific Use:

Knitted cord used to retract the gingival area before placing impression material

Product Name:

Gingival retraction cord non-impregnated dental sz0 10' cotton

Special Features:

Size 0; used w/astringent liquid, pn 0111; recommended for anterior teeth; used as bottom cord where two cord technique is desired

Material:

Cotton

Material:

Cotton

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T243-e