

View Online at <https://aerobasegroup.com/nsn/6520-01-272-6068>

**Length:**

10.000 feet

**Physical Form:**

Cord

**Size Designation:**

2

**Disposition After Initial Use:**

Disposable

**Specific Use:**

General purpose knitted cord used to retract the gingival area before placing impression material

**Product Name:**

Gingival retraction cord non-impregnated used as gp cord

**Special Features:**

Used w/astringent liquid, pn 0111; recommended for bicuspid and molars; used as an upper cord for those who desire a two-cord technique

**Material:**

Cotton

**Material:**

Cotton

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T243-e