

View Online at <https://aerobasegroup.com/nsn/6520-01-335-4648>

**Compression Strength:**

8500.0 pounds per square inch

**Usage Design:**

Denture purpose and model

**Setting Time:**

Between 3 minutes and 5 minutes

**Weight:**

50.0 pounds

**Product Name:**

Stone artificial dental powder form pail 50 lb

**Special Features:**

Powder form; working time 10 to 60 seconds; setting expansion 0.08 pct; mixing ratio 26 to 28 ml water to 100 gm powder; used to secure denture models to articulator mounting plates; for model stone when a rapid set is required; for obtaining interocclusal records and for locking various parts of gnathological instruments (articulators, pantographs etc) into position; unit qty pail 50 lb.

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T243-d