

View Online at <https://aerobasegroup.com/nsn/6520-01-335-4648>

Compression Strength:

8500.0 pounds per square inch

Usage Design:

Denture purpose and model

Setting Time:

Between 3 minutes and 5 minutes

Weight:

50.0 pounds

Product Name:

Stone artificial dental powder form pail 50 lb

Special Features:

Powder form; working time 10 to 60 seconds; setting expansion 0.08 pct; mixing ratio 26 to 28 ml water to 100 gm powder; used to secure denture models to articulator mounting plates; for model stone when a rapid set is required; for obtaining interocclusal records and for locking various parts of gnathological instruments (articulators, pantographs etc) into position; unit qty pail 50 lb.

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d