

View Online at <https://aerobasegroup.com/nsn/6520-01-366-7546>

**Physical Form:**

Powder

**Usage Design:**

Denture restoration

**Weight:**

500.0 grams

**Special Features:**

Methylmethacrylate copolymer powder; cold cure method; for mixing w/acrylic liquid; mixing ratio: 1 0 g powder to 7 ml liquid; handling and curing time subject to room temp. Of 23 degrees c; used for finishing of chrome cobalt partials and functional repairs such as reconstruction of margins and base lining

**Color:**

Pink, light red

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T243-d