

View Online at <https://aerobasegroup.com/nsn/6520-01-366-7546>

Physical Form:

Powder

Usage Design:

Denture restoration

Weight:

500.0 grams

Special Features:

Methylmethacrylate copolymer powder; cold cure method; for mixing w/acrylic liquid; mixing ratio: 1 0 g powder to 7 ml liquid; handling and curing time subject to room temp. Of 23 degrees c; used for finishing of chrome cobalt partials and functional repairs such as reconstruction of margins and base lining

Color:

Pink, light red

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d