

View Online at <https://aerobasegroup.com/nsn/6520-01-464-8482>

Physical Form:

Liquid and powder

Usage Design:

Denture restoration

Weight:

21.0 grams

Curing Process Type:

Heat

Special Features:

Characterized lucitone denture base resin; 250 powders, 21 grams each; 10 gingival tone powders, 10 grams each; with accessories; yields 250 units

Color:

Pink, light

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d