

View Online at <https://aerobasegroup.com/nsn/6520-01-481-9124>

Physical Form:

Liquid

Usage Design:

Denture restoration

Setting Time:

4 minutes

Curing Process Type:

Self

Product Name:

Resin restorative dental liquid denture restoration 4 minutes

Special Features:

Temp-bond clear, temporary resin cement, for use in cementation of temporary and provisional restoration such as crowns, bridges, inlays and outlays; flexible and transparent when cured; includes a 3 gram tube of base, a 3 gram tube of catalyst, and a mixing pad

Color:

Clear

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T243-d