

View Online at <https://aerobasegroup.com/nsn/6525-01-064-7334>

Cooling Method:

Oil

Rectification Type:

Full-wave

Tube Insert Anode Type:

Stationary

Tube Insert Radiographic Peak Voltage Rating:

100.0 kilovolts

Err Default:

Err default

Err Default:

Err default

Focal Spot Quantity:

1

Focal Spot Size:

1 mm square

Transformer Rectification Method:

Solid state

Err Default:

Err default

Product Name:

Tube-transformer head dental x-ray appar stationary 100 kilovolt

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

T140-c