

View Online at <https://aerobasegroup.com/nsn/6525-01-070-1547>

**Mounting Type:**

Mobile

**Accessory Components And Quantity:**

Installation, calibration, maintenance, operator manual sets 2; field use bucky table 1

**Current Rating In Amps:**

75.0

**Design Type:**

Mobile

**Operation Method:**

Manual

**Control Type:**

Cabinet

**Tube Quantity:**

1

**Electric Power Rating:**

25000.0 watts

**Tube Type:**

Double focus, rotating anode radiographic

**Table Type:**

Nontilting

**Rectification Type:**

Full-wave

**Tube Cooling Method:**

Air

**Tube Mounting Type:**

Tube stand

**Tube Stand Mounting Type:**

Spring counter-balance

**Radiographic Timer Type:**

Exposure

**Transformer Rectification Method:**

Solid state

**Transformer Output Capacity:**

250 ma at 100 kv peak and 200 ma at 125 kv peak

**Product Name:**

X-ray app rad med mobile 250ma at 100pkv 115-230 volt 50-60hz ac

**Voltage In Volts:**

115.0 and 230.0

**Phase:**

Single

**Frequency In Hertz:**

50.0 or 60.0

**Supplementary Features:**

C/o collimeter, control unit, table, transformer, and tube stand

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

Yes - demil/mli

**Fig:**

T140-c