

View Online at <https://aerobasegroup.com/nsn/6525-01-456-9987>

**Physical Form:**

Liquid

**Unit Package Yield:**

7.000 liters

**Material Processing Method:**

Automatic

**Processing Speed:**

Rapid

**Unit Package Content Form:**

Single unit for mixing

**Special Features:**

Developer-replenisher, 6 liter bottle, yields to 7 liters, 2 bottles per case; special package chemicals designed for 40 seconds processing in the fuji medical laser imager film special containers of concentrate to make 7 liters of working solution are designed to fit the belt end chemical mixer in the film integrated processor

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A24700