

View Online at <https://aerobasegroup.com/nsn/6525-01-456-9987>

Physical Form:

Liquid

Unit Package Yield:

7.000 liters

Material Processing Method:

Automatic

Processing Speed:

Rapid

Unit Package Content Form:

Single unit for mixing

Special Features:

Developer-replenisher, 6 liter bottle, yields to 7 liters, 2 bottles per case; special package chemicals designed for 40 seconds processing in the fuji medical laser imager film special containers of concentrate to make 7 liters of working solution are designed to fit the belt end chemical mixer in the film integrated processor

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A24700