X-ray Film Processing Fixer - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6525-01-531-1045

Physical Form:   Liquid   Quantity Within Each Unit Package:   4.240 liters   Material Processing Method:   Manual   Unit Package Content Form:   Premixed unit   Product Name:   Fixer x-ray film processing fixer pg 4.24 li   Special Features:   Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray   processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case   Shelf Life:   N/a   Unit Of Measure:
Quantity Within Each Unit Package:         4.240 liters         Material Processing Method:         Manual         Unit Package Content Form:         Premixed unit         Product Name:         Fixer x-ray film processing fixer pg 4.24 li         Special Features:         Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case         Shelf Life:         N/a         Unit Of Measure:
4.240 liters Material Processing Method: Manual Unit Package Content Form: Premixed unit Product Name: Fixer x-ray film processing fixer pg 4.24 li Special Features: Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
Material Processing Method: Manual Unit Package Content Form: Premixed unit Product Name: Fixer x-ray film processing fixer pg 4.24 li Special Features: Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
Manual Unit Package Content Form: Premixed unit Product Name: Fixer x-ray film processing fixer pg 4.24 li Special Features: Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
Unit Package Content Form: Premixed unit Product Name: Fixer x-ray film processing fixer pg 4.24 li Special Features: Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
Premixed unit Product Name: Fixer x-ray film processing fixer pg 4.24 li Special Features: Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
Product Name: Fixer x-ray film processing fixer pg 4.24 li Special Features: Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
Fixer x-ray film processing fixer pg 4.24 li Special Features: Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
Special Features: Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
Insta-neg manual x-ray processing is reduced to less than one minute; premixed for use in either standard developing tanks for total x-ray processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
processing including panoramic film or in small containerts for endodontic and emergency processing; (4) developers and (4) mixers per case Shelf Life: N/a Unit Of Measure:
case Shelf Life: N/a Unit Of Measure:
Shelf Life: N/a Unit Of Measure:
N/a Unit Of Measure:
Unit Of Measure:
Demilitarization:
No
Fiig:

A24700