## NSN 6530-00-079-7039

Sterilization Indicator - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6530-00-079-7039

Physical Form:
Strip
Product Name:
Indicator sterilization indicator color chg by steam 250s
Special Features:
Indicator color change by steam heat; center line changes color from white to brown when proper sterilization conditions of steam;
temperature and time have occured within pack, at 250 degf aafter 12 to 15 min. Exposure, and 270 degrees f after 2 to 3 min exposure;
steam; item w/potency period or expiration date; unit qty, 250 strips
Shelf Life:
N/a
Unit Of Measure:
Demilitarization:
No
Fiig:
T396-a