## NSN 6530-00-181-7770

Sterilization Indicator - Page 1 of 1

T396-a



View Online at https://aerobasegroup.com/nsn/6530-00-181-7770 **Physical Form:** Sheet **Unit Package Type:** Box **Product Name:** Indicator sterilz sensitized paper or fabric test record9x12"31s **Exterior Container Weight:** 33.00 pounds **Exterior Container Cubic Measure:** 1.641 cubic feet **Exterior Container Quantity:** 48 **Special Features:** Flexible to conform to irregular contours in surgical packs, a light background shall be over-printed w/a darker colored steam-heat sensitive ink, the ink shall undergo a uniform color change in the presence of a vacuum, or absence of air when exposed to steam heat for 3-1/2 min. At 272 degrees f, or for 15 min at 250 degrees f, unit q type 31 **Intermediate Package Quantity:** 12 **Unit Of Issue Cube:** 0.010 cubic feet **Unit Of Issue Weight:** 0.35 pounds Shelf Life: N/a **Unit Of Measure: Demilitarization:** No Fiig: