NSN 6530-01-163-7302

Sterilization Indicator - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6530-01-163-7302

| _ | | | |
|----|-------|------|--------|
| Sn | acial | Feat | IIPAC' |
| | | | |

| Sensitized paper; | color changes by steam heat; | white bar will turn black | when proper | sterilization ha | as been achieved; | 2 minute exposure |
|--------------------|-------------------------------------|----------------------------|-----------------|------------------|-------------------|-------------------|
| time at 270 deg. F | ; 5/8 in. By 8 in. (long form); for | or determining the efficie | ncy of steriliz | ation process | | |

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T396-a