NSN 6540-01-302-1024

Ophthalmic Perimeter - Page 1 of 1

T146-j



View Online at https://aerobasegroup.com/nsn/6540-01-302-1024

| Current Type: |
|--|
| Ac |
| Accessory Components And Quantity: |
| Bowl, hemispherical-1; crt display monitor w/electric light pen mode selector-1; built-in printer-1; color filters- (red, blue, green) |
| Design Designation: |
| Goldman |
| Product Name: |
| Perimeter ophthalmic goldman ac 110/220v 50/60hz |
| Special Features: |
| 110/220 volt; 50/60 hz; stimulus intensity range equal to or greater than 0.08 to 10, 000 asb; 33 cm bowl radius w/bowl illumination of 31.5 |
| asb; custom design also stores color patterns; vertically adjustable; statistical software pkg for comparing 10 previous exams and |
| calculates significant changes; 14 cubic ft volume, under 150 lb; used for performing test for glaucoma screening, central 30 screening, |
| peripheral screening, central threshold and neurological threshold |
| Shelf Life: |
| N/a |
| Unit Of Measure: |
| - |
| Demilitarization: |
| No |
| Fiia· |